



A decision with a future

DRÄGER EVITA®
INFINITY® V500

How do I ensure staying prepared for the future?



Decision for high quality ventilation therapy

There is no singular method to fulfilling one of the most critical challenges in today's ICU. Reduction of the patient's time on a ventilator remains a key issue. A high quality level of ventilation therapy at all times needs to be maintained through a consistent standard of treatment plans.



Everyday you are faced with new patients, new challenges. Minutes count when assessing a patient's condition to formulate a diagnosis and take the right decision on therapy.

Once the patient is stabilized, when is the right moment to start weaning? A clinician knows that weaning from a mechanical ventilator can take an extended period of time. How can this be done as quickly as possible?

THE EVITA STANDARD

Evita devices have supported you for over 25 years with

- high quality standards
- configurability and upgradeability
- advanced training and service concepts

Meeting your daily challenges now and in the future will require modern-day solutions. The Evita V500 ventilator offers innovative technology and is synonymous with providing high performance ventilation.

A complete set of treatment tools

TOOLS FOR YOUR VENTILATION THERAPY

- Lung diagnostic tools like the Low Flow maneuver (inflection points)
- Breath by breath recruitment trends (e.g. PEEP, EIP, VT, Cdyn)
- Recruitment tools (e.g. Inspiration Hold, QuickSet, PressureLink)
- PC-APRV with AutoRelease
- Volumetric CO₂-Monitoring (VCO₂, VTCO₂, Slope Phase 3, Vds/VTe)
- Weaning parameter (e.g. RSBi, P0.1, NIF)
- Automatic weaning with SmartCare/PS®
- Variable Pressure Support, Proportional Pressure Support, Automatic Tube Compensation
- Graphical representation of airway resistance and lung compliance with Smart Pulmonary View
- "Room-to-breathe" concept (AutoFlow®, BIPAP, VG)
- Applicable for neonatal ventilation (smallest tidal volume: 2 ml)


FUNCTIONS TO SUPPORT YOUR WORKFLOW

- Non-invasive ventilation in all modes and all patient categories
- O₂-Therapy allows flow constant oxygen application
- Flexible screen configuration: 6 different views for each patient related to the individual therapy
- Full record of all patient data, alarms and trends
- Data export via USB interface
- Context-sensitive help function and online Instruction For Use for ventilation modes, alarms and device functions
- RFID functionality e.g. to monitor exchange intervals and to transport ventilation settings



D-76108-2013

Sophisticated support in every therapy phase.



The Dräger Evita V500 accompanies you along the respiration pathway and delivers high level therapy and decision support – for your daily challenges today and tomorrow.

PREVENT

- Non-invasive ventilation with sensitive inspiratory triggering and expiratory cycling-off to improve patient-ventilator synchrony

STABILIZE

- Low Flow analysis of static lung characteristics to optimize the ventilation settings
- QuickSet and PressureLink allow to facilitate virtually any recruitment procedure
- AutoRelease in PC-APRV stabilizes the end-expiratory lung volume

WEAN

- SmartCare/PS as an automated, knowledge-based weaning software reduces length of stay in the ICU*
- Configurable limits in SmartCare/PS allow individual adaptations to patient situation
- Comprehensive monitoring function to observe the progress of weaning
- Variable Pressure Support to increase the natural tidal volume variability

RECOVER

- O₂-Therapy supports recovery after extubation and non-invasive ventilation

* F. Lellouche et al.; Am J respir Care Med Vol 174, pp 894-900, 2006. Results are based on a European Multicenter Randomized Trial with 144 patients demonstrating improved respiratory condition, with stable hemodynamic and neurologic status, and no ARDS prior to initiating weaning

Transport



D-43483-2012

GS500, PS500, TRANSPORT SUPPLY UNIT (TSU), BED DOCKING

Flexible use almost anywhere in the hospital for up to four hours of independent ventilation thanks to the external gas and power supplies. The transport supply unit for two gas cylinders and the bed coupling make the patient transport easy.

Accessories



MT-2099-2007

INFINITY ID ACCESSORIES

The Infinity® ID accessory components simplify use and optimize workflows: The RFID functions detect faulty connections and coupled systems, and can import ventilation settings when connected to the ventilator. Replacement intervals of accessories can be monitored.

Training



D-76030-2013

EVITA V500 PRODUCT TRAINER

Dräger supports your daily hospital work with innovative training programs. Our interactive product trainers provide useful information about our products as well as practical instruction about the intuitive operating system and use of the device.

Service



D-32436-2011

DRÄGERSERVICE – BECAUSE QUALITY COUNTS

For any requirement, we have the service contract that fits: from inspection-only and preventive maintenance to complete service packages. We have the answer for your individual needs.

For further information and downloads please visit www.draeger.com

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany

www.draeger.com

Manufacturer:

Dräger Medical GmbH
Moislinger Allee 53–55
23558 Lübeck, Germany

REGION EUROPE CENTRAL AND EUROPE NORTH

Dräger Medical GmbH
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION EUROPE SOUTH

Dräger Médical S.A.S.
Parc de Haute Technologie d'Antony 2
25, rue Georges Besse
92182 Antony Cedex, France
Tel +33 1 46 11 56 00
Fax +33 1 40 96 97 20
dlmfr-contact@draeger.com

REGION MIDDLE EAST, AFRICA

Dräger Medical GmbH
Branch Office
P.O. Box 505108
Dubai, United Arab Emirates
Tel +971 4 4294 600
Fax +971 4 4294 699
contactuae@draeger.com

REGION ASIA / PACIFIC

Draeger Medical
South East Asia Pte Ltd.
25 International Business Park
#04-27/29 German Centre
Singapore 609916, Singapore
Tel +65 6572 4388
Fax +65 6572 4399
asia.pacific@draeger.com

REGION NORTH AMERICA

Draeger Medical, Inc.
3135 Quarry Road
Telford, PA 18969-1042, USA
Tel +1 215 721 5400
Toll-free+1 800 437 2437
Fax +1 215 723 5935
info.usa@draeger.com

REGION CENTRAL AND SOUTH AMERICA

Dräger Panama Comercial S. de R.L.
Complejo Business Park,
V tower, 10th floor
Panama City
Tel +507 377 9100
Fax +507 377 9130
contactcsa@draeger.com

As of August 2015:

Dräger Medical GmbH changes
to Drägerwerk AG & Co. KGaA.