

KIMAL VASCULAR ACCESS

INTERNATIONAL
PRODUCT
PORTFOLIO

PATIENT SAFETY

CRITICAL CARE

PAGE 4

ICU
Operating Theatre
HDU
Emergency
Cardiac



AT THE
HEART

IV ACCESS

PAGE 8

Oncology
Haematology
OPAT
Radiology
Nutrition

RENAL

PAGE 10

Haemodialysis
Peritoneal Dialysis
Outpatient (Satellite Centres)
Radiology

OF WHAT WE DO


ACCESSORIES

PAGE 12

Renal
Critical Care
IV Access


ALTIUS® CVC

- LuerSafe™ automatic clamps to reduce risk of air emboli and infection.
- Colour-coded lumens to reduce misidentification of lines.
- Including high pressure lumen for pressures of up to 300psi.
- Superior insertion kit for safe and smooth insertion.



ALTIUS PROACTIV+ CVC

- Protecting your patients for the life of the catheter*
- ProActiv+ non-eluting antimicrobial
- Chlorhexidine and silver sulfadite-free.
- No leaching chemicals released into the bloodstream
- Superior insertion kit for safe and smooth insertion.



ALTIUS SEVEN CVC

- The only 7 lumen CVC at 8.5 French.
- Propriety Multi-tube™ technology provides laminar flow for superior flow rates.
- Superior insertion kit for safe and smooth insertion.



ALTIUS RT CATHETER

- Peace of mind from air embolism, severe blood loss and catheter kinking.
- LuerSafe automatic clamps.
- Soft, flexible PUR material.
- Excellent flow rates at low pressures.
- Superior insertion kit for safe and smooth insertion.



KFLOW EPIC ST DIALYSIS CATHETER

- Excellent flow rates at low pressures.
- Thermoresponsive radiopaque polyurethane.
- Modified Double D design for excellent performance.
- Superior insertion kit for safe and smooth insertion.



KIMAL TRANSDUCERS

- Venous and arterial transducers available as venous only, arterial only, and arterial and venous combined.
- Deltran Utah disposable transducer.
- Compatible with most monitor devices.



KIMAL MIDLINE CATHETER

- Protect your patients' long-term vessel health.
- Clearly identified MIDLINE to prevent accidents.
- For therapies of more than 10 to 30 days.
- Two different insertion kits to meet the needs of different inserters.
- Available in 1 to 3 lumens, and 4 to 6 French.



KIMAL PICCLINE CATHETER

- Safe, reliable long-term access, available in 1-3 lumens.
- CT-rated up to 300psi 5ml/sec, identified with a dedicated purple lumen.
- Contains two nitinol guidewires (45cm and 80cm) for clinician preference.
- Two different insertion kits to meet the needs of different inserters.



ARTERIAL PACKS

Available with customised insertion kits and arterial catheters pack.



VASPACKS INSERTION PACKS

- Create your own maximum barrier for catheter insertion.
- Put your Kimal catheter inside your pack.
- Bespoke arterial pack with pink buttoned tourniquets.
- High quality components and no unexpected surprises.
- Fast turn around from as little as 6 weeks from design to delivery.



ALTIUS GUIDEWIRE

- Superior nitinol guidewire with depth and ECG markings.
- Unique ergonomic dispenser.
- Thumb-controlled dispenser, even with soiled gloves.
- From paediatric to adults



ACCESSORIES

- We offer a comprehensive range of accessories for your infusion and therapies; from manometers, to clamps and spikes. Check our accessories section to review our most popular accessories.



IV ACCESS

NUTRITION

RADIOLOGY

OPAT

HAEMATOLOGY

ONCOLOGY

KIMAL PICCLINE CATHETER

- Safe, reliable long-term access, available in 1-3 lumens.
- CT-rated up to 300psi 5ml/sec, identified with a dedicated purple lumen.
- Contains two nitinol guidewires (45cm and 80cm) for clinician preference.
- Two different insertion kits to meet the needs of different inserters.



KIMAL MIDLINE CATHETER

- Protect your patients' long-term vessel health.
- Clearly identified MIDLINE to prevent accidents.
- For therapies of more than 10 to 30 days.
- Insertion kit to meet the needs of different inserters.
- Available in 1 to 3 lumens, and 4 to 6 French.
- Two different insertion kits to meet the needs of different inserters.



IV ACCESS

ONCOLOGY

HAEMATOLOGY


OPAT

RADIOLOGY

NUTRITION

KIMAL MST KITS

- Create your own unique insertion kit with nitinol guidewires, sharps safety components and but-toned tourniquets.



VASPACKS INSERTION PACKS

- Create your own maximum barrier for catheter insertion.
- Put your Kimal catheter inside your pack.
- Bespoke arterial pack with pink buttoned tourniquets.
- High quality components and no unexpected surprises.
- Fast turn around from as little as 6 weeks from design to delivery.



ALTIUS RT CATHETER

- Peace of mind from air embolism, severe blood loss and catheter kinking.
- LuerSafe automatic clamps.
- Soft, flexible PUR material.
- Excellent flow rates at low pressures.
- Superior insertion kit for safe and smooth insertion.



KFLOW EPIC ST DIALYSIS CATHETER

- Excellent flow rates at low pressures.
- Thermoresponsive radiopaque polyurethane.
- Modified Double D design for excellent performance.
- Superior insertion kit for safe and smooth insertion.



KFLOW EPIC LT DIALYSIS CATHETER

- Split and stepped tip designs
- Thermoresponsive radiopaque polyurethane
- Superior insertion kit for safe and smooth insertion.
- Adult and paediatric range available.



KFLOW PD CATHETER

- Reliable PD catheter allowing you to build your own unique insertion pack.



VASPACKS INSERTION PACKS

- Create your own maximum barrier for catheter insertion.
- Put your Kimal catheter inside your pack.
- Bespoke arterial pack with pink buttoned tourniquets.
- High quality components and no unexpected surprises.
- Fast turn around from as little as 6 weeks from design to delivery.



KIMAL ON/OFF & FISTULA PACKS

- Customised packs to support renal therapy.
- High quality components and no unexpected surprises.
- Fast turn around from as little as 6 weeks from design to delivery.



K39574 On/off pack


K39575 Fistula pack



**Low pressure infusion line
700cm**
K39/10/700
Box qty: 50



Manometer line plain
K39/39 | K39/40 | K39/41 |
K39/42 | K39/120
Available in 50, 100, 150 and
200cm length.
Box qty: 50



Manometer line colour strip
Available in red or blue colour
strip, and in 100, 150 and
200cm.
Red: K39/30, K39/31 and
K39/32
Blue: K39/35, K39/36 and
K39/37
Box qty: 100



**Badge fistula
compression bandages**
DR8/2
Box qty: 200



Connector
K402
Box qty: 1,400



**Female to female
connector**
K408
Box qty: 150



Perforator
K405
Box qty: 100



Priming set
K39/8
Box qty: 100



Clamp scissors blue
K61/10
Box qty: 200




Peritoneal catheter adaptor
K48/2
Box qty: 100




Trocar
KSTROCAR
Box qty: 1



**Connector valve | Repair
valve**
KFE/COV | KFE REV
Box qty: 5



Y-Site accessory set
K46125
Box qty: 30



Y-site accessory set
K407
Box qty: 350



Transducer protector sterile
K34/3/LL
Box qty: 300



**Pressure monitoring kit
single**
K200/1
Box qty: 25



**Double pressure transducer
set**
K200/4
Box qty: 20




Arterial Transducer Set
K200/1/RE
Box qty: 25

NOTES

NOTES

**WHY
CHOOSE
KIMAL?**

SUPPLYING HOSPITALS
IN MORE THAN 
70 countries
worldwide

**OVER
50
YEARS**
experience

**TREATING
+2M**
patients
EVERY YEAR



All delivered with the
**KIMAL
COMMITMENT**

DELIVER ON TIME | HIGH QUALITY PRODUCTS | NO UNEXPECTED SURPRISES
AWARD-WINNING CUSTOMER SERVICE | SPECIALIST CLINICAL KNOWLEDGE

For more information, pricing or samples, speak to your local Kimal Representative, or contact:

INTERNATIONAL

Tel: +44 (0)845 437 9543
Fax: +44 (0)845 437 9544
Email: internationalenquiries@kimal.co.uk

MIDDLE EAST + AFRICA

Tel: +97 16 552 8146
Fax: +97 16 552 8147
Email: mea@kimal.co.uk

© 2019 Kimal PLC. All rights reserved. No part of this document may be reproduced in any form or by any means without permission in writing from Kimal PLC. Kimal PLC reserve the right to change the specification of its products at any time. Kimal PLC, Arundel Road, Uxbridge, UB8 2SA.

DCR 1458 A1 | JULY 2019

KIMAL
delivering healthcare innovation