



# Mini Care

## Compressor Nebulizer

Medication delivery can be a pleasant experience! Mini Care, a new breed of APEX Nebulizer, is designed to make medication administration more palatable to the daily life and further encourage patients to comply with the prescribed treatments. This piston-driven unit has been designed to deliver swift and tranquil treatments. Its low noise level also makes it practical for late-night treatment. Mini Care is an ultra-portable unit with a compact design. It is equipped with a prominent carrying handle, a recessed tubing connection for durability and an adequate storage space for easy access of tubing and other accessories. Most importantly, Mini Care assures that patients' medication will be delivered safely, efficiently, and reliably every time at any time.

Care for a Healthy Life

**APEX**

# Mini Care

## Compressor Nebulizer

### Features ▼

- All-in-one packing (everything packs up into main unit).
- Convenient built-in handle.
- Large storage compartment (for the nebulizer kit).
- Hidden electrical cord compartment.
- Convenient tube holder (for the nebulizer kit).
- Appealing exterior design.
- Reliable, proven performance - CE/FDA 510(k) certified.



Large compartment area to store tubes and accessories.



Power cord can be neatly tucked into the back of the unit.



Built-in nebulizer holder offers hands free medication delivery.



**Mouthpiece set :** Anti-drool "T" mouthpiece with 6" reservoir tube and 7' supply tube.

**Mouthpiece set and 2 masks :** Anti-drool "T" mouthpiece with 6" reservoir tube, 7' supply tube and one of each adult & pediatric mask

### Specification

Electrical Rating	120 VAC, 60 Hz, 1.2 A
(Note: Refer to the rating label on the product)	230 VAC, 60 Hz, 0.6 A
	230 VAC, 50 Hz, 0.6 A
Maximum Compressor Pressure	≥ 30 psi
Operation Flow Rate	≥ 5.0 lpm
Nebulization Rate	0.2 - 0.3 ml / min
Particle size	0.5 - 5.0 microns
Dimensions (W x D x H)	29.8 × 17.3 × 10.7 cm / 11.7" × 6.8" × 4.2"
Weight	Approx. 1.9 kg / 4.2 lb

Distributed by

Apex Medical Corp.  
 No. 9, Min Sheng St., Tu-Cheng,  
 New Taipei City, 23679, Taiwan  
 T +886 2 2268 3100  
 F +886 2 2268 6525  
 marketing@apexmedicalcorp.com  
 www.apexmedicalcorp.com  
 Print-2014/07 / All rights reserved  
 Specifications may change without notice.

